

UltraFreezeTM

Liquid nitrogen cryosurgical system

The compact, lightweight
UltraFreeze is the new,
easy-to-use, portable system
ideal for most office
skin procedures.

The newly designed ergonomic
handle offers comfort and control.


Completely self-contained, the
system delivers controlled spray.

UltraFreeze is supplied with an assortment
of easily attached luer lock closed-tip
or open apertures.

健力醫療器材有限公司
CLINICON MEDICAL LTD.
Unit 411, Lippo Sun Plaza,
28 Canton Road,
Tsimshatsui, Kowloon,
Hong Kong

Tel: 2730 1883

Fax: 2736 1106



Nothing works like a **WALACH**

A cost effective cryosurgical system with Wallach's proven reliability

Totally enclosed to
protect user from cold

Luer-Lock connector for
quick change of apertures

New smooth flow
spray technology

Rotating trigger for
ambidextrous use

Closed tips
for small lesions

Ergonomically designed
grip for comfort and control

**The new Wallach UltraFreeze is the ideal
choice for a wide variety of common disorders.**

Common Verrucae

Periungual Warts

Plantar Warts

Actinic Keratoses

Seborrheic Keratoses

Lentigines


Small Keloids or Hypertrophic Scars

Leukoplakia

Sebaceous Hyperplasia

Cavernous Hemangiomas

Available in 0.5
and 0.3 liters



AccuShield

*Plastic reusable plate with
6 conical openings. Helps to direct liquid
nitrogen spray in a concentrated area.*

New plated brass apertures
available in a variety of sizes

*For additional information call 1-800-243-2463
or visit our website: www.wallachsurgical.com*

Nothing works like a **WALLACH**