

10

Series 2 35 Watt

Flexible/gooseneck arm

L 220 212 A without stand

L 220 222 A ¹⁾ on stand with caster-base

L 220 232 A ¹⁾ on height adjustable stand with caster-base

Colour of lamp white

Bulb type Halogen with cool-light reflector

Bulb details 12 V / 35 W, 10°

Light intensity appr. 22.500 Lux / 500 mm

Light field Ø 90 / 180 mm
at 500 / 1.000 mm

Colour temperature appr. 3.000 (3.700)² Kelvin

Life time appr. 3.000 h (bulb)

Switch at the lamp head

Working radius appr. 800 mm

Electrical connection 230 V, 50 - 60 Hz

Total weight appr. 3,2 kg / 9,1 kg¹

Joint pin 18 x 35 mm

Stand height appr. 1.070 mm

Caster-base steel, white coated

Base diameter appr. 635 mm

Casters Twin 50 mm plastic casters,
(2 x antistatic and lockable, 3 x normal)

optional: Dimmer / international plug /

² colour temperature filter / 110 V - transformer



doubled-wall grounding warmth
directing construction

high value aluminium cool- body

