



DR-008FD Operating Stool



Specifications
Lifting Stroke 195 mm
Net Weight 17 kg
Upholstery Choices --- 025

- Seat height adjustment by gas cylinder.
- Smooth and light seat height adjustment by ring lever.
- Highly stabilized and compact base for operating rooms and treatment rooms.
- Foot control keeps hands clean.


DR-007FD Operating Stool



Specifications
Lifting Stroke 195 mm
Net Weight 18 kg
Upholstery Choices --- 025

- Seat height adjustment by gas cylinder.
- Smooth and light seat height adjustment by ring lever.
- Highly stabilized and compact base for operating rooms and treatment rooms.
- 360° rotatable backrest and 250mm adjustable height for the ideal posture of the doctor.
- Foot control keeps hands clean.


DR-010DN



Specifications
Lifting Stroke 180 mm
Net Weight 29 kg
Rating Voltage
AC 110/ 220/ 230/ 240 V
Upholstery Choices --- 025

- Seat height adjustment by electro-hydraulic pump.
- Smooth and light seat height adjustment by pedal switch.
- Low initial height and height adjustment for the ideal posture of the doctor.


HY-1200 Operating Stool



Specifications
Lifting Stroke 200 mm
Net Weight 27 kg
Upholstery Choices --- M1

- Seat height adjustment by hydraulic pump.
- Height adjustable armrest with sliding movement.
- Foot control keeps hands clean.
- Stable and compact base for easy access to surgical areas.

HY-1100 Operating Stool



Specifications
Lifting Stroke 200 mm
Net Weight 18 kg
Upholstery Choices --- 025
Backrest Option

- Seat height adjustment by hydraulic pump.
- Foot control keeps hands clean.
- Stable and compact base for easy access to surgical areas.

The manufacturer reserves the right to alter specifications without prior notice.

Belmont®

TAKARA BELMONT CORPORATION
2-1-1, Higashishinsaibashi,
Chuo-ku, Osaka, 542-0083, Japan
TEL : 81-6-6213-5945
FAX : 81-6-6212-3680